

Using the *DEBUSSY v2.0beta* with the Graphical User Interface: a quick guide

Contents

Introduction.....	2
Quick guide for the <i>CLAUDE_GUI</i>	2
<i>Building the sampled interatomic distances of a population of NCs</i>	3
<i>Building the sampled interatomic distances for a single cluster</i>	4
<i>Building the sampled interatomic distances starting from an .xyz file</i>	5
<i>Computing a powder pattern from a single cluster (.smp file)</i>	6
<i>Computing the Pair Distribution Function from a powder diffraction pattern</i>	7
Quick guide for the <i>DEBUSSY_GUI</i>	8

Introduction

The *DebUsSy Suite* is an open source software package implementing the Debye scattering equation (DSE) for the analysis of powder diffraction pattern from nanocrystalline, defective and/or non-periodic materials.

The *DebUsSy Suite* strategy is of bottom-up type and consists of two main steps, managed respectively by the *CLAUDE* Suite (generating populations of NCs) and the Debussy program (implementing the calculation of DSE and the model parameters optimization against the experimental data).

This quick guide provides a general overview of the multiple tasks that can be performed using the *DebUsSy Suite*, release v2.0 (both *CLAUDE* and Debussy programs), through the Graphical User Interface (*GUI*). The users can directly perform some of these tasks. A **tutorial** is available at <http://debussy.sourceforge.net/debussy/tutorial.html>, including input files, experimental data and detailed instructions.

Quick guide for the *CLAUDE_GUI*


The *CLAUDE_GUI* can perform several tasks, which are shortly introduced here, in the following order:

- Building the sampled interatomic distances of a population of NCs.
- Building the sampled interatomic distances for a single cluster.
- Building the sampled interatomic distances starting from a user-defined .xyz file.
- Computing a powder pattern from a single cluster (.smp file).
- Computing the Pair Distribution Function from a powder diffraction pattern.

All the *CLAUDE_GUI* tasks require selecting the *CLAUDE* tab of the *GUI* (1).

All the operations described in each section must be performed in the order indicated by the numbered lists.

Building the sampled interatomic distances of a population of NCs



2) Select your working folder by clicking on **Set Folder** button

3) Select the **POPULATION DATABASE** tab.

4) Upload your **.cif/.pha** file by clicking on **LOAD .cif/.pha** button and select the right file. The **.cif** file is automatically converted in the **.pha** file that is displayed in the **GUI** window. Edit this **.pha** file for re-organizing the number/type of atomic species, according to the specificity of your structure, if necessary. Otherwise all the atoms with the same atomic symbol will be grouped in a single atomic species.

5) Fill the **POPULATION DATABASE** tab or load an already existent **.ddb** file using **LOAD .ddb** button.

6) **SAVE** the **.ddb** file.


7) Click on **CELL** button of the **MAKE** vertical menu (the unit cell is built starting from the asymmetric unit and from the space group information reported in the **.pha** file).

8) Click on **DATABASE** button of the **MAKE** vertical menu. The **DISTANCES/SAMPTO*** folder is created in the working folder, including the **.smp** and **.smp_INFO** files of each cluster.

9) Optionally, click on the **XYZ** button of the **MAKE** vertical menu to generate the **.xyz** files containing the atomic Cartesian coordinates of the NCs. It is strongly recommended to limit the maximum number of the NCs for this operation, by editing the **.ddb** file.

10) If step 9 has been performed, click on the **XYZ** button of the **PLOT** vertical menu in order to visualize and manage the atomistic models of the NCs (**Jmol** application will be launched).

Building the sampled interatomic distances for a single cluster




Follow the list of operations from 1 to 5 reported for building the population of NCs.

6) Select the **largest_only** option in the **POPULATION DATABASE** tab.

7) **SAVE** the **.ddb** file.


8) Click on **LARGEST** in the **MAKE** vertical menu to generate a single cluster database (single **.smp/smp_INFO** in the case of a monovariate population, single base size combined with all the possible lengths along **c** for bivariate populations of NCs), having the largest size specified in the **.ddb** file.

Building the sampled interatomic distances starting from an .xyz file



- 2) Select your working folder by clicking on **Set Folder** button.
- 3) Select the **SINGLE CLUSTER DATABASE** tab:
- 4) Upload the .xyz file of the cluster by clicking on the **...** button.
- 5) Fill the **SINGLE CLUSTER DATABASE** tab with the information about the phase and the sampling (see the *Claude User's Manual* for details) or load an existing .ddb file by clicking on **LOAD .ddb**.
- 6) **SAVE** the .ddb file.
- 7) Click on the **MOLEC** button of the **MAKE** vertical menu. The DISTANCES/SAMPTO* folder is created in the working folder, including a .smp and a .smp_INFO files.

Computing a powder pattern from a single cluster (.smp file)



2) Select your working folder by clicking on **Set Folder** button.

3) Select the **PATTERN** tab.


4) Fill the **PATTERN** tab or load an existing *diffractor.inp* file by clicking on **LOAD .inp** button and edit it, if necessary.

5) **SAVE** the *diffractor.inp* file.

6) Click on the **PATTERN** button of the **MAKE** vertical menu, in order to compute the pattern from the .smp indicated in the *diffractor.inp* file.

7) The output file with the .tqi extension is created and can be visualized automatically using the **PATTERN** button of the **PLOT** vertical menu.

Computing the Pair Distribution Function from a powder diffraction pattern




- 2) Select your working folder by clicking on **Set Folder** button.
- 3) Select the **PDF** tab.
- 4) Fill the **PDF** tab or load an existing *dopdf.inp* file by clicking on **LOAD .inp** button and edit it, if necessary.
- 5) **SAVE** the *dopdf.inp* file.
- 6) Click on the **PDF** button of the **MAKE** vertical menu. Several files can be created, according to the options selected in the input file (see the *Claude User's Manual* for details)
- 7) In order to automatically visualize the reduced $G(r)$ computed, the **G(r)** button in the **PLOT** vertical menu can be used.

Quick guide for the *DEBUSSY_GUI*

The Debussy program can be used for the following tasks:

- simulating a powder diffraction pattern by the DSE in the absence or in the presence of an experimental pattern;
- refining a model powder pattern against the experimental data;
- calculating the e.s.d. of the refined parameters.

A database of sampled interatomic distances is necessary for any of the previous listed tasks (see the *CLAUDE_GUI* section).



- 1) Select the *DEBUSSY* tab of the *DebUsSy Suite GUI*.
- 2) Select your working folder by clicking on **Set Folder** button.
- 3) Fill the required input files (*.dwa*, *.par*, *.ref*) using the related tabs or load the existing ones using the **LOAD .dwa/.par/.ref** buttons.
- 4) Click on the **Simulation** button of the **RUN** vertical menu in order to compute the model pattern described by the set of parameters written in the *.par* file(s). This action does not perform any refinement of parameters and can be used to preliminary testing the starting model.
- 5) Click on the **Refinement** button of the **RUN** vertical menu in order to perform the refinement of the model parameters.
- 6) Click on the **STD** button of the **RUN** vertical menu, compute the e.s.d. of the refined parameters, once the minimum of the χ^2 function is reached.
- 7) Use the **PLOT** buttons in the vertical menu to visualize all the graphical outputs generated during and at the end of the refinement or the simulation:

Simulation to plot the pattern calculated by the **Simulation** run.

Best fit to plot the pattern calculated by the **Refinement** run

Size to plot the mono- or bivariate size distributions.

Lattice exp. to plot the lattice parameter(s) as a function of the NCs size.

S.O.F. to plot the atomic site occupancy factors as a function of the NCs size.

BTH to plot the atomic Debye Waller factors as a function of the NCs size.